



EMG/NCV/EP  
Model: 5000Q  
PRO Portable  
**Negar Andishgan**

## EMG/NCV/EP Model: 5000Q PRO Portable 4 Channels

**Accurate  
Reliable  
Easy to use**

Negar Andishgan Manufacturing **EMG/NCV/EP5000Q Portable** is a real portable device, full featured diagnostic system, a low noise high performance amplifier and state of the art software features.

- Portable EMG system with Light weigh and easy to port
- 2 channel Low noise amplifier
- Designed to achieve the high flexibility, portability.
- Easy to setup and State of the Art software Development.
- Generating customize report in world, PDF, or Video.
- Compact and efficient user interface
- Attractive price





**Accurate  
Reliable  
Easy to use**

**Model:5000Q  
EMG/NCV/EP**





**Protocol wizard:**  
This is a powerful tool for defining new tests, designing and customizing your system in your own way.

Name of nerves, muscles, number of signal trace, method of calculation, shape of result Table, final report and many other options can be defined by Protocol wizard.



- Microsoft® Word report generator can help you to generate your report in less than a minute.
- Graphical directory for muscle and nerve electrode placement can give you useful recommendations.
- There are many kind of applicable functions such as:
  - Left and Right comparison
  - Super imposer
  - Pick detector
  - Needle Video and Audio reviewer
  - Signal annotation to help you get more precise results.



- Microsoft® SQL Server Guarantees the integrity of data management which can be accessed trough network.



Address: No.19, Koshesh Ave, North Sohrevardi, Tehran ,Iran  
Tel: +982188529469 - 70 Fax: +982188529563